

OHICHO III

Balloon Compliance (Typical value)

Balloon Diameter (mm)	Pressure (atm [$\times 10^5$ Pa])											
	2	4	6	8	10	12	14	16	18	20	22	24
4.0	3.55	3.63	3.73	3.83	3.92	4.00	4.07	4.11	4.16	4.20	4.24	<i>4.31</i>
5.0	4.43	4.51	4.65	4.78	4.89	5.00	5.07	5.14	5.19	5.25	5.31	<i>5.39</i>
6.0	5.25	5.42	5.60	5.77	5.89	6.00	6.08	6.15	6.21	6.28	<i>6.36</i>	<i>6.48</i>
7.0	6.03	6.25	6.49	6.70	6.86	7.00	7.12	7.22	7.31	<i>7.44</i>	<i>7.58</i>	—

 Nominal Pressure

 Rated Burst Pressure

※Do not apply pressure exceeding the Rated Burst Pressure.

※Italic font is used for the reference values.