

Balloon Compliance (Typical value)

Balloon Diameter (mm)	Pressure (atm [$\times 10^5$ Pa])										
	2	4	6	8	10	12	14	16	18	20	22
1.50	1.33	1.37	1.40	1.44	1.47	1.50	1.52	1.55	1.58	1.60	
1.75	1.55	1.60	1.63	1.68	1.72	1.75	1.77	1.81	1.84	1.87	
2.00	1.72	1.78	1.84	1.90	1.96	2.00	2.03	2.05	2.08	2.11	
2.25	1.94	1.99	2.07	2.14	2.20	2.25	2.29	2.33	2.35	2.38	2.42
2.50	2.21	2.27	2.35	2.40	2.46	2.50	2.54	2.58	2.62	2.66	2.71
2.75	2.37	2.43	2.51	2.60	2.68	2.75	2.82	2.87	2.91	2.94	2.97
3.00	2.63	2.70	2.77	2.86	2.93	3.00	3.05	3.09	3.13	3.16	3.19
3.25	2.87	2.96	3.05	3.13	3.20	3.25	3.30	3.34	3.37	3.40	3.45
3.50	3.08	3.17	3.26	3.36	3.44	3.50	3.56	3.61	3.65	3.68	3.72
3.75	3.29	3.39	3.49	3.59	3.68	3.75	3.81	3.86	3.90	3.94	3.99
4.00	3.50	3.61	3.72	3.83	3.92	4.00	4.06	4.11	4.15	4.20	4.26
4.25	3.67	3.83	3.97	4.09	4.17	4.25	4.31	4.36	4.41	4.47	
4.50	3.86	4.04	4.20	4.31	4.41	4.50	4.56	4.62	4.67	4.73	

Nominal Pressure
 Rated Burst Pressure

※Do not apply pressure exceeding the Rated Burst Pressure.

※Italic font is used for the reference values.