

Balloon Compliance

Balloon Diameter (mm)	Pressure atm/kPa																					
	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27
	608	709	811	912	1013	1115	1216	1317	1419	1520	1621	1723	1824	1925	2027	2128	2229	2330	2432	2533	2634	2736
2.00	1.84	1.87	1.90	1.93	1.96	1.98	2.00	2.01	2.03	2.04	2.05	2.06	2.08	2.10	2.11	2.13	2.15	2.17	2.19	2.22	-	-
2.25	2.07	2.10	2.14	2.17	2.20	2.23	2.25	2.27	2.29	2.31	2.33	2.34	2.35	2.37	2.38	2.40	2.42	2.44	2.46	2.48	2.50	2.53
2.50	2.35	2.38	2.40	2.43	2.46	2.48	2.50	2.52	2.54	2.56	2.58	2.60	2.62	2.64	2.66	2.69	2.71	2.74	2.77	2.79	2.82	2.85
2.75	2.51	2.56	2.60	2.64	2.68	2.72	2.75	2.79	2.82	2.85	2.87	2.89	2.91	2.93	2.94	2.95	2.97	2.99	3.00	3.03	3.06	3.09
3.00	2.77	2.81	2.86	2.90	2.93	2.96	3.00	3.03	3.05	3.07	3.09	3.11	3.13	3.15	3.16	3.17	3.19	3.20	3.22	3.26	3.29	3.32
3.25	3.05	3.09	3.13	3.17	3.20	3.23	3.25	3.28	3.30	3.32	3.34	3.35	3.37	3.38	3.40	3.43	3.45	3.48	3.50	3.54	3.57	3.61
3.50	3.26	3.32	3.36	3.40	3.44	3.47	3.50	3.53	3.56	3.58	3.61	3.63	3.65	3.67	3.68	3.70	3.72	3.74	3.76	3.81	3.85	3.89
3.75	3.49	3.55	3.59	3.64	3.68	3.72	3.75	3.78	3.81	3.83	3.86	3.86	3.90	3.92	3.94	3.96	3.99	4.02	4.04	4.09	4.14	4.18
4.00	3.72	3.78	3.83	3.88	3.92	3.96	4.00	4.03	4.06	4.08	4.11	4.11	4.14	4.18	4.20	4.23	4.26	4.30	4.33	4.37	4.42	4.46
4.25	3.97	4.03	4.09	4.13	4.17	4.21	4.25	4.28	4.31	4.33	4.36	4.39	4.41	4.44	4.47	4.50	4.53	4.57	4.62	4.66	-	-
4.50	4.20	4.26	4.31	4.37	4.41	4.46	4.50	4.53	4.56	4.59	4.62	4.46	4.67	4.70	4.73	4.76	4.79	4.84	4.89	4.94	-	-

 Nominal Pressure(NP)

 Rated Burst Pressure(RBP)

*Do not apply a pressure exceeding RBP value

(Note)Italics:reference value