

Compliance Chart

| Balloon Diameter (mm) | Pressure (atm) | | | | | | | | | | | | | | | | | | |
|-----------------------|----------------|------|------|------|------|------|------|------|------|------|------|------|------|-------------|-------------|-------------|-------------|-------------|-------------|
| | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 |
| 2.00 | 1.76 | 1.80 | 1.85 | 1.89 | 1.94 | 1.97 | 2.00 | 2.03 | 2.06 | 2.08 | 2.10 | 2.12 | 2.15 | <i>2.17</i> | <i>2.19</i> | <i>2.21</i> | <i>2.23</i> | <i>2.25</i> | <i>2.25</i> |
| 2.50 | 2.17 | 2.27 | 2.33 | 2.39 | 2.44 | 2.47 | 2.50 | 2.53 | 2.55 | 2.57 | 2.59 | 2.62 | 2.64 | <i>2.66</i> | <i>2.68</i> | <i>2.71</i> | <i>2.74</i> | <i>2.77</i> | <i>2.79</i> |
| 3.00 | 2.57 | 2.69 | 2.79 | 2.84 | 2.89 | 2.96 | 3.00 | 3.04 | 3.07 | 3.10 | 3.13 | 3.16 | 3.19 | <i>3.22</i> | <i>3.25</i> | <i>3.28</i> | <i>3.31</i> | <i>3.35</i> | <i>3.39</i> |
| 3.50 | 3.02 | 3.16 | 3.25 | 3.33 | 3.39 | 3.45 | 3.50 | 3.54 | 3.57 | 3.60 | 3.64 | 3.67 | 3.70 | <i>3.73</i> | <i>3.77</i> | <i>3.80</i> | <i>3.84</i> | <i>3.88</i> | <i>3.92</i> |
| 4.00 | 3.53 | 3.67 | 3.75 | 3.83 | 3.89 | 3.95 | 4.00 | 4.04 | 4.07 | 4.11 | 4.15 | 4.19 | 4.22 | <i>4.26</i> | <i>4.29</i> | <i>4.34</i> | <i>4.39</i> | <i>4.43</i> | <i>4.48</i> |

Nominal Pressure

Rated Burst Pressure: **Do not apply a pressure exceeding this value.**

(Note) Italics; reference value