

# Compliance Chart

Balloon Size (mm)	Pressure (atm)																						
	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
2.50	2.18	2.22	2.25	2.29	2.32	2.36	2.40	2.43	2.45	2.48	2.50	2.51	2.53	2.55	2.56	2.57	2.58	2.59	2.61	2.62	2.64	<i>2.65</i>	<i>2.67</i>
2.75	2.36	2.41	2.45	2.49	2.54	2.58	2.62	2.66	2.68	2.72	2.75	2.77	2.79	2.81	2.83	2.85	2.86	2.87	2.88	2.89	2.91	<i>2.93</i>	<i>2.95</i>
3.00	2.60	2.67	2.71	2.75	2.80	2.84	2.88	2.91	2.94	2.97	3.00	3.02	3.04	3.06	3.08	3.09	3.11	3.12	3.14	3.15	3.17	<i>3.19</i>	<i>3.20</i>
3.25	2.82	2.89	2.94	2.99	3.04	3.08	3.13	3.16	3.19	3.22	3.25	3.28	3.30	3.32	3.33	3.35	3.36	3.38	3.40	3.42	3.44	<i>3.46</i>	<i>3.48</i>
3.50	3.03	3.10	3.16	3.22	3.27	3.32	3.37	3.41	3.43	3.47	3.50	3.53	3.55	3.57	3.58	3.60	3.62	3.63	3.65	3.67	3.70	<i>3.72</i>	<i>3.75</i>
3.75	3.25	3.33	3.39	3.45	3.51	3.56	3.61	3.65	3.68	3.72	3.75	3.78	3.80	3.82	3.84	3.87	3.88	3.90	3.92	3.94	3.97	<i>4.00</i>	<i>4.03</i>
4.00	3.47	3.56	3.63	3.69	3.75	3.80	3.86	3.90	3.93	3.97	4.00	4.02	4.05	4.07	4.10	4.13	4.15	4.17	4.19	4.21	4.23	<i>4.25</i>	<i>4.28</i>
4.25	3.67	3.76	3.83	3.91	3.97	4.03	4.09	4.13	4.17	4.21	4.25	4.28	4.31	4.33	4.36	4.39	4.41	4.44	4.47	<i>4.50</i>	<i>4.53</i>	-	-
4.50	3.86	3.97	4.04	4.12	4.20	4.26	4.31	4.37	4.41	4.46	4.50	4.53	4.56	4.59	4.62	4.64	4.67	4.70	4.73	<i>4.76</i>	<i>4.79</i>	-	-

  Nominal Pressure

  Rated Burst Pressure: **Do not apply a pressure exceeding this value.** (Note) Italics; reference value