

Balloon Compliance (Typical value)

Balloon Diameter (mm)	Pressure (atm [$\times 10^5$ Pa])								
	2	4	6	8	10	12	14	16	18
4.0	3.62	3.77	3.89	4.00	4.07	4.16	4.23	<i>4.31</i>	<i>4.41</i>
5.0	4.60	4.79	4.92	5.00	5.06	5.11	5.17	<i>5.25</i>	<i>5.37</i>
6.0	5.21	5.46	5.75	6.00	6.15	6.26	6.37	<i>6.48</i>	<i>6.61</i>

 Nominal Pressure

 Rated Burst Pressure

※Do not apply pressure exceeding the Rated Burst Pressure.

※Italic font is used for the reference values.