

Compliance Chart

Balloon Size (mm)	Pressure (atm)									
	4	5	6	7	9	11	13	15	16	17
3.00	2.84	2.91	2.97	3.00	3.10	3.17	3.25	3.32	<i>3.36</i>	<i>3.39</i>
4.00	3.69	3.77	3.85	4.00	4.02	4.10	4.18	4.26	<i>4.29</i>	<i>4.35</i>
5.00	4.76	4.86	4.95	5.00	5.16	5.29	5.43	5.59	<i>5.67</i>	<i>5.77</i>
6.00	5.74	5.87	5.99	6.00	6.26	6.43	6.62	6.83	<i>6.96</i>	<i>7.09</i>
7.00	6.60	6.76	6.90	7.00	7.21	7.39	7.58	7.79	<i>7.91</i>	<i>8.03</i>
8.00	7.36	7.56	7.74	8.00	8.11	8.32	8.53	8.74	<i>8.86</i>	<i>9.00</i>
9.00	8.12	8.36	8.60	9.00	9.05	9.29	9.52	<i>9.78</i>	<i>9.94</i>	<i>10.11</i>
10.00	9.15	9.42	9.67	10.00	10.17	10.46	10.76	<i>11.10</i>	<i>11.25</i>	<i>11.42</i>

 規定圧力 (Nominal Pressure)
 定格破裂圧 (Rated Burst Pressure)

※定格破裂圧を超えて使用しないでください。 *斜体字: 参考値を示す