

取引先各位

製造販売元：株式会社カネカ  
販売元：株式会社カネカメディックス

**カネカエクスチェンジデバイス CO-N1 『KUSABI』  
8Fr ガイディングカテーテル適合のご案内**



拝啓 時下ますますご清栄のことと、お慶び申し上げます。  
平素は弊社製品に格別のご高配を賜り、厚く御礼申し上げます。  
さて、弊社製品カネカエクスチェンジデバイス CO-N1 『KUSABI』につきまして、発売より 6・7Fr ガイディングカテーテル適合としておりましたが、本年 8 月より 8Fr ガイディングカテーテルにつきましても追加適合となりましたのでご案内申し上げます。なお製品仕様に変更が無いため、今後は 6・7・8Fr ガイディングカテーテル用としてご使用いただけます。また拡張圧、カタログ番号、JAN コードに変更は御座いません。今後とも、より一層のご愛顧を賜りますよう、宜しくお願い申し上げます。

敬具

記

- 製品名 : KUSABI (クサビ)
- 販売名 : カネカエクスチェンジデバイス CO-N1
- 承認番号 : 22400BZX00362000
- 一般名称 : 冠動脈カテーテル交換用カテーテル (肺動脈カテーテル交換用カテーテル)
- 対象製品 : カタログ番号 KS-107、KS-107LB、KS-117
- 適用時期 : 2020年8月より

- 主な変更内容 : ①ガイディングカテーテル適合サイズに 8Fr が追加。  
②製品ラベル内の推奨ガイディングカテーテルサイズに 8Fr を追加 (赤枠)。

＜変更前＞	＜変更後＞
 <p>The image shows the product label for KUSABI before the change. The 'Guidance Catheter Compatibility' section lists '6Fr or 7Fr'. The 'Catheter Effective Length' is 1070 mm. The 'Rated Burst Pressure' is 14 atm. The label includes various technical specifications and safety instructions.</p>	 <p>The image shows the product label for KUSABI after the change. The 'Guidance Catheter Compatibility' section now lists '6Fr 7Fr 8Fr', with '8Fr' highlighted in a red box. The 'Catheter Effective Length' is 1070 mm. The 'Rated Burst Pressure' is 14 atm. The label includes various technical specifications and safety instructions.</p>

以上