

Compliance Chart

Balloon Diameter (mm)	Pressure (atm)																													
	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	
4.0	3.92	3.97	4.00	4.03	4.07	4.08	4.08	4.08	4.09	4.09	4.09	4.10	4.10	4.10	4.11	4.11	4.11	4.11	4.11	4.11	4.11	4.11	4.11	4.11	4.12	4.12	4.12	4.13	4.13	
5.0	4.89	4.95	5.00	5.04	5.07	5.08	5.09	5.09	5.10	5.10	5.10	5.11	5.11	5.11	5.12	5.12	5.13	5.13	5.13	5.13	5.14	5.14	5.14	5.15	5.15	5.15	5.16	5.16	5.17	
6.0	5.86	5.92	6.00	6.05	6.07	6.07	6.09	6.10	6.11	6.11	6.11	6.12	6.12	6.12	6.13	6.13	6.14	6.14	6.15	6.16	6.16	6.17	6.17	6.18	6.18	6.19	6.19	6.20	6.20	
7.0	6.81	6.90	7.00	7.06	7.09	7.10	7.11	7.11	7.12	7.13	7.13	7.14	7.15	7.16	7.17	7.17	7.18	7.18	7.19	7.19	7.20	7.20	7.21	7.21	7.22	7.22	7.23	7.23	7.23	

Nominal Pressure
 Rated Burst Pressure

※Do not apply pressure exceeding the Rated Burst Pressure.