

# Compliance Chart

Balloon Size (mm)	Pressure (atm)																	
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
1.50	-	1.37	1.39	1.41	1.43	1.46	1.49	1.50	1.52	1.53	1.55	1.56	1.57	1.59	<i>1.60</i>	<i>1.61</i>	<i>1.62</i>	<i>1.64</i>
2.00	-	1.76	1.81	1.86	1.91	1.94	1.98	2.00	2.03	2.05	2.07	2.09	2.11	2.13	<i>2.15</i>	<i>2.17</i>	<i>2.18</i>	<i>2.21</i>
2.25	-	1.97	2.03	2.09	2.14	2.18	2.21	2.25	2.28	2.30	2.32	2.35	2.37	2.40	<i>2.43</i>	<i>2.45</i>	<i>2.48</i>	<i>2.51</i>
2.50	-	2.25	2.30	2.35	2.40	2.44	2.47	2.50	2.53	2.55	2.58	2.60	2.63	2.65	<i>2.68</i>	<i>2.70</i>	<i>2.73</i>	<i>2.77</i>
2.75	-	2.49	2.55	2.60	2.65	2.69	2.72	2.75	2.77	2.80	2.82	2.85	2.88	2.90	<i>2.93</i>	<i>2.96</i>	<i>2.99</i>	<i>3.02</i>
3.00	-	2.66	2.74	2.81	2.87	2.93	2.96	3.00	3.03	3.06	3.10	3.12	3.15	3.18	<i>3.21</i>	<i>3.25</i>	<i>3.29</i>	<i>3.32</i>
3.25	-	2.91	2.99	3.06	3.12	3.17	3.21	3.25	3.28	3.31	3.35	3.37	3.40	3.44	<i>3.47</i>	<i>3.50</i>	<i>3.53</i>	<i>3.57</i>
3.50	-	3.13	3.21	3.29	3.35	3.41	3.45	3.50	3.53	3.57	3.60	3.64	3.68	3.72	<i>3.78</i>	<i>3.83</i>	<i>3.86</i>	<i>3.91</i>

Nominal Pressure

Rated Burst Pressure: **Do not apply a pressure exceeding this value.**

(Note) Italics; reference value