


Balloon Compliance (Typical value)

Balloon Diameter (mm)	Pressure (atm [$\times 10^5$ Pa])								
	2	4	6	8	10	12	14	16	18
1.00	0.97	0.98	1.00	1.02	1.05	1.07	1.09	<i>1.11</i>	<i>1.14</i>
1.50	1.44	1.47	1.50	1.53	1.55	1.58	1.61	<i>1.65</i>	<i>1.69</i>
2.00	1.86	1.93	2.00	2.05	2.09	2.12	2.16	<i>2.21</i>	<i>2.25</i>
2.25	2.08	2.16	2.25	2.31	2.35	2.40	2.44	<i>2.48</i>	<i>2.53</i>
2.50	2.29	2.39	2.50	2.57	2.62	2.67	2.72	<i>2.76</i>	<i>2.81</i>
2.75	2.53	2.64	2.75	2.83	2.88	2.93	2.99	<i>3.04</i>	<i>3.10</i>
3.00	2.77	2.88	3.00	3.09	3.15	3.20	3.25	<i>3.32</i>	<i>3.38</i>

 Nominal Pressure

 Rated Burst Pressure

- ※Do not apply pressure exceeding the Rated Burst Pressure.
- ※Italic font is used for the reference values.