

Compliance Chart

Balloon Diameter (mm)	Pressure (atm)												
	2	3	4	5	6	7	8	9	10	11	12	13	14
2.50	2.35	2.38	2.42	2.48	2.50	2.54	2.56	2.59	2.60	2.63	2.65	2.67	2.69
3.00	2.81	2.86	2.91	2.96	3.00	3.04	3.07	3.11	3.13	3.16	3.18	3.21	3.24
3.50	3.26	3.31	3.37	3.44	3.50	3.55	3.60	3.64	3.67	3.70	3.73	3.77	3.79
4.00	3.72	3.75	3.85	3.93	4.00	4.06	4.12	4.17	4.21	4.25	4.28	4.32	4.35

Nominal Pressure
 Rated Burst Pressure

※Do not apply pressure exceeding the Rated Burst Pressure.

※Italic font is used for the reference values.