

Balloon Compliance

Ballon Diameter (mm)	Pressure (atm/kPa)														
	2	4	6	8	10	12	14	16	18	20	22	24	26	28	30
	203	405	608	811	1013	1216	1419	1621	1824	2027	2229	2432	2634	2837	3040
2.0	1.79	1.80	1.82	1.84	1.86	1.89	1.91	1.94	1.97	2.00	2.03	2.05	2.08	2.11	2.15
2.5	2.18	2.19	2.21	2.23	2.25	2.28	2.33	2.38	2.44	2.50	2.55	2.60	2.66	2.71	2.77
3.0	2.56	2.58	2.61	2.64	2.66	2.74	2.80	2.87	2.94	3.00	3.06	3.12	3.18	3.25	3.32
4.0	3.28	3.31	3.35	3.40	3.47	3.55	3.66	3.78	3.90	4.00	4.10	4.19	4.29	4.40	4.50
5.0	4.00	4.03	4.09	4.18	4.30	4.46	4.62	4.76	4.88	5.00	5.12	5.26	5.40	5.56	5.73
6.0	4.61	4.68	4.81	4.99	5.20	5.38	5.54	5.68	5.83	6.00	6.17	6.37	-	-	-

Normal Pressure (NP) Rated Burst Pressure (RBP) (Note) Italics:reference value

* A balloon diameter of 6.0 mm indicates the NP and RBP are identical

Do not apply a pressure exceeding RBP are identical